



QuikPath RAPID

Point-of-care Isothermal PCR for Covid-19

10

viral copies/ μ L

Limit of Detection

95%

Sensitivity



2 Gene Detection
with internal control




HSA Approved



Safe
Zero Aerosol



Scalable
Scale-up without
additional capex




Adaptable
Diverse use cases



Traceable
Barcode detection

Product Specifications

“Plug and Play” cartridges combined with basic heat block allow for affordable scaling up of POC testing



- 20°C storage



Results in 40 minutes




Result interpretation with colorimetry



Uses basic heat block



6 months shelf life



Contact Us

sg.lifescience@biomed-global.com
www.quikpath.sg